

Battery Eliminator



MAIN FEATURES:

Battery Eliminator is designed to be safe and reliable source of low voltage, DC power. It is ideal for use in classrooms and substitutes for a wide variety of batteries including AA, AAA, C, D, and 9-volt types. Allows teachers to conduct experiments requiring a variety of direct voltages while avoiding the decreased characteristics of dry cells. Comprised of a heavy duty transformer wound with copper wire with high lamination to avoid losses. Voltage is easily selected with the turn of a dial. Multiple tapping outputs of 2, 4, 6, 8, 10, 12 Volts DC.

SPECIFICATIONS:

- Product Code : GE 104 2-12 V DC/1 Amp.
- Product Code : GE 105 2-12 V DC/2 Amp.
- Product Code : GE 106 2-12 VDC/3 Amp.
- Product Code : GE 107 2-12 V DC/4 Amp.
- Product Code : GE 108 2-12 V DC/5 Amp.