

IC Regulated Battery Eliminators (Stabilized)



MAIN FEATURES:

This is a reliable and safe source of low D.C. voltage. This regulated Battery Eliminator is regulated through solid state LC.'s which makes its output highly stable and hence perfect for laboratory use. The design of the circuit and the construction ensures that the Battery Eliminator continues to deliver accurate & stable performance over its long life.

SPECIFICATIONS:

- Product Code : GE 110 2-12 V/1 A in the steps of 2 Volts.
- Product Code : GE 111 2-12 V/2 A in the steps of 2 Volts.
- Product Code : GE 112 2-12 V/3 A in the steps of 2 Volts.
- Product Code : GE 113 2-12 V/5 A in the steps of 2 Volts.