

Photo diode Characteristics Apparatus

Product Code:GE 232

Objective: To Plot a Graph, Current Vs intensity of Light at different value of Voltages



Specifications:

- Instrument comprises of DC Regulated Power Supply 0-3 VDC/150mA,
- Round Meters for Voltage & Current Measurement,
- Photo diode mounted in a Metal box,
- Connections of Supplies & Photo Diode brought out at 4mm Sockets.
- Lamp holder with 60W Bulb & 50 cm Wooden Plank are given with this apparatus.