

LDR Characteristics Apparatus

Product Code :GE 230

Objective: To Plot a Graph, Current Vs Intensity of Light at Different Value of Voltages.



Specifications:

- Instrument comprises of DC Regulated Power Supply 0-3 VDC/150mA,
- Two Round Meters For Voltage & Current Measurement,
- LDR mounted in a Metal Box,
- Connections of Supplies & LDR brought out at 4mm Sockets.
- Lamp Holder with 60 W Bulb & 50 cm Wooden PLank are given with this apparatus