

## Digital Flux Meter / Gauss Meter

Productcode: GE490



The Gaussmeter operates on the principle of Hall Effect in Semiconductors. The small Hall voltage is amplified through a high stability amplifier that a millivoltmeter connected at the output of the amplifier can be calibrated directly in magnetic Field unit.

### Technical Specifications:

- Range: 0-2K gauss and 0-20 Kg
- \* Resolution: One Gauss at 0-2 Kgauss range
- Accuracy :  $\pm 0.5\%$
- Display: 3 1/2 digit? seg. (DPM)
- Detector; Hall probe with an imported hall element
- Power; 220V. 50 Hz
- Special: Indicate the direction of the magnetic field.