

## Magnetic Pickup Trainer

Product code: GE 805



### Specifications:

- \* Study of Speed Measurement by variable reluctance Transducer
- 12V, 2400 RPM permanent magnet DC motor
- Speed controller (Q-2400 RPM)
- Electromagnetic pickup
- Signal conditioner circuitry (OPAMP based)
- Excitation voltage For pickup
- \* 4 digit digital counter for speed display
- Test points at various places
- Electronic tacho generation
- Interconnection: 2/4mm banana sockets
- DC Supply: Built in IC regulated power supplies
- 220V  $\pm 10\%$ , 50Hz mains operated
- Enclosed in an attractive ABS plastic cabinet
- Accessories Users with patch cords