

## Angular Displacement Trainer

Product Code: GE810



### Specifications:

- Servo potentiometer as angular displacement measurement transducer (full 0-3600) using calibrated dial with 2' resolution
- Linear potentiometer as linear displacement measurement transducer (0-100mm) using mechanical attachment with vernier
- Opamp in signal conditioner circuit
- Offset adder for position modulation
- 3.5 digit DPM to read displacement
- Constant calibrated
- Interconnection: 2/4 mm banana sockets
- DC Supply: Built in IC regulated power supply
- 220V  $\pm 10\%$ , 50Hz mains operated
- Enclosed in an attractive ABS plastic cabinet with cover
- Accessories: Users Manual with patch cords