


Applications OF IC 555

Product Code: GE341

Objective : Kit for IC 555 Timer circuit applications has been designed to perform the following experiments:-

- Astable Multivibrator
- Monostable Multivibrator
- Bistable Multivibrator
- Voltage to Time convertor / Voltage to Frequency convertor



PTechnical Specifications:

- Inbuilt Fixed DC Regulated Power supply Output Voltages: +5VDC
- IC & Components Provided IC: 555
- Potentiometer: 1 No.
- Switches (Push to ON): 2 Nos.
- Different resistance & capacitors are provided on front panel.
- Power requirement: 230 VAC 10%. 50Hz.
- Power chord. Patch chords & Instruction manual.