


## LVDT Transducer trainer

Product code: GE 802



### Specifications:

- Measurement Range 10/20 mm
- Excitation Frequency: 2,5 KHz (approximately)
- Sensitivity: 10mV DC/mm
- Linear Range: Full Scale
- Signal conditioner output: 0.1 VDC or Maximum Displacement
- Display: 3 1/2% Digit LED with Polarity Indicator
- Micrometer Least count: 0.01 mm
- inbuilt IC regulated DC supply
- Power Supply: 220 V,  $\pm 1$  0%, 50 Hz
- Standard Accessories: Users Manual with patch cords