


Measurement of Temperature Using RTD Trainer (PT 100)

Product Code :GE 803



Specifications:

- Study of RTD as Temperature measuring Transducer
- Mercury thermometer and heater 10W (optional)
- RTD(PT-100) probe with protect cover
- Bridge circuit For transducer
- Differential amplifier with feedback
- Lead compensation system
- Sockets provided to study the process to find out the error rate
- 3.5 digit DPM as temperature indicator
- Interconnection: 2/4mm banana sockets
- 220V \pm 10%, 50HZ mains operated
- Enclosed in an attractive ABS plastic cabinet with cover
- Accessories: User's Manual with patch cords