

## DC Speed Control System Trainer

Product Code: GE 913



### Technical Specifications

- Study of speed control of DC motor in open and close loops.
- Speed control of 12V, 4W, permanent mag net DC motor.
- Speed range upto 3000 RPM (typical).
- Electronic techo- generator for feedback.
- Non-contact eddy current break for load.
- 4 digit RPM counter.
- 3YdigitDVM.
- Built in regulated power supply.
- 220V  $\pm 10\%$ , 50Hz mains operated.
- User's Manual with patch cords