

IV-E228D Michealson Interferometer

List of Items Supplied With Setup:

Sr.No.	Experimental set up Contain	Quantity
1	Optical Breadboard (450 X 600 mm)	1
2	Kinematic laser mount	1
3	Mirror mount with translation	1
4	Mirror mount with precision translation	1
5	Rotation stage (Resolution 10/division)	1
6	Screen with mount (Dimension 75x75mm)	1
7	Pressure cell (Pressure Range 0-300 mm Hg)	1
8	Mirror	2
9	Beam splitter mount	1
10	Beam splitter	1
11	Glass slide	5
12	Beam Expander with mount	1
13	Diode laser 5mW (650nm or 532nm)	1
14	Diode laser power supply	1
15	Post Base	4