

IV-E262 Figure Of Merit Of A Ballistic Galvanometer Trainer



Scope of Learning:

- DETERMINATION OF FIGURE OF MERIT OF A BALLISTIC GALVANOMETER

Technical Specification:

- Power Supply: +2V DC, at 100mA
- Variable Resistor Box: 1-1K, 1-100K, 1-100K (P, Q and R)
- Push Key: One
- 4Way Key: One
- Mains: 230V AC $\pm 10\%$, 50Hz
- Fuse: 500Ma
- Dimension of Unit (mm): W 298 x D 248 x H 130

Salient Feature:

- Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.
- Ballistic Galvanometer With Lamp and Scale Arrangement.
- Instruction manual.
- Connections are brought out through 4mm Colored Sockets.
- Patch Cords 4mm.
- The trainer is housed in ABS Plastic cabinet.
- Size of the trainer set 12"x10"

Optional Accessories:

- No