

## IV-E337A BH Curve Tracer Trainer



### Scope of Learning:

- Hysteresis loss

### Specification:

- Three Different Values of Capacitors are given on board.
- One Variable Potentiometer
- Different Circuit diagram is printed on board.
- Step Down Transformer with Demonstration transformer with 2 secondary coils and one Primary Coil.
- Output through BNC Connectors and BNC Lead are Provided for Output to CRO
- Components mounted on the front panel.

### Features:

- Trainer operated on 230V, 50 Hz + 10%
- Supplied with instruction manual
- Supplied with sufficient nos. of patch cords
- Supplied with Computer Lead
- Dimensions of the Trainer will be 12"X8" with Plastic Box Highly shock proof.