

## IV-E383 Digital Gauss Meter



- Gauss Meter is an instrument to find the Magnetic Field of Any Object

### Experiment:

- To find the Magnetic Flux (Electromagnetic Field of any Object)

### Digital Gauss Meter:

- Range: 0-2 K Gauss & 0-20 K Gauss,
- Resolution: 1 Gauss at 0-2 K Gauss Range,
- Accuracy:  $\pm 0.5\%$ ,
- Display:  $3\frac{1}{2}$  Digit, 7 Segment LED,
- Power: 220V, 50Hz
- Special Feature: Indicates the direction of the Magnetic field.
- Gauss Probe

### Salient Features:

- Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.
- Instruction manual.
- Connections are brought out through 3Pin Connector
- The trainer is housed in ABS Plastic and Metal Cabinet.
- Size of the trainer set 12"x4.5"

### Optional Accessories:

- Electromagnet.