



GURU KRIPA ENTERPRISES

AN ISO 9001:2015 & NSIC CERTIFIED COMPANY

Sur No - 668, Plot No - 268 , Fatevadi, Sarkhej , Ahmedabad,Gujarat - 382210

Email: info@gkelabs.com,www.gkelabs.com,Mobile::9106316361



GEM PORTAL
REGISTERED

IV-E388C Photo-Electric Effect Current VS. Intensity & Wavelength Of Light; Max. Energy Of Photo-Electrons VS Frequency Of Light

Salient Feature:

- Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.
- Instruction manual.
- Connections are brought out through BNC Connector.
- The trainer is housed in ABS Plastic cabinet/Metal Box.
- Size of the trainer set 12"x5"

Packing List:

- Wooden Optical Bench With Scale Stickered
- Set of 5 Color Filter
- Photo Cell Unit Housed in Plastic Box
- Incandescent Lamp Housed in Aluminum Box
- BNC to BNC Probe
- Photo Electric Effect Trainer Unit
- Lux Meter (Optional)

