



# GURU KRIPA ENTERPRISES

AN ISO 9001:2015 & NSIC CERTIFIED COMPANY

Sur No - 668, Plot No - 268 , Fatevadi, Sarkhej , Ahmedabad,Gujarat - 382210

Email: info@gkelabs.com,www.gkelabs.com,Mobile::9106316361



GEM PORTAL  
REGISTERED

## IV-E592D To Determine The Wavelength Of He-Ne Laser Using Diffraction Grating

### Apparatus Supply:

- He-Ne Laser Wavelength: 630nm, Output 1mW
- Diffraction Grating: 500L/inch
- Single Slit
- Double Slit
- Single Wire
- Cross Wire
- Mesh Wire and Some More
- Analyzer Polarizer
- Quarter Wave Plate
- Half Wave Plate
- Photo Detector With Photo Sensor
- Optical Bench Rail Type/Rod Type 1.5mth.
- Optical White Screen

