


## IV-E723 Young Modulus Of A Wire By Optical Lever Method



### Scope of Learning:

- YOUNG MODULUS OF A WIRE BY OPTICAL LEVER METHOD

### Required Accessories:

- A Long assembly with optical lever mirror, sample wire.
- Red Light Pointer.
- Inch Tape.
- Telescope with scale arrangement
- Frame holding the experimental wire
- Clamping cylinder with hanging arrangement
- Set of slotted weights
- Meter Rod
- Screw Gauge

### Optional Accessories:

- No