

IV-E968 Constant Current Power Supply



Constant Current Power Supply:

(i) Digital Mill Ammeter

- Range: 0-1A
- Accuracy: $\pm 0.1\%$ of reading ± 1 digit

(ii) Constant Current Power Supply:

- Current 0-1A
- Resolution: $100\mu\text{A}$
- Accuracy: $\pm 0.2\%$ of the reading ± 1 digit
- Load regulation: 0.03% for 0 to full load
- Line regulation: 0.05% for 10% variation

(iii) Termination:

- 2Nos. (4mm for Current)

Salient Features:

- Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.
- Instruction manual.
- Connections are brought out through 3Pin Connector
- The trainer is housed in ABS Plastic and Metal Cabinet.
- Size of the trainer set 12"x4".