

## IV-MS100S Magnetic Stirrer

- Digital LCD display. Control through Knob.
- Maximum Stirring Volume 2L.
- Stirring - MS100S.
- Magnetic stirrer is mainly suitable for liquid mixing, it can simultaneously do heating and stirring of low viscosity liquids or solid liquid mixtures. Its used as an experimental tool in various Universities, Environmental Protection, Sanitation, Epidemic Prevention, Chemical Tap Water, Medical and Other Units.



IV-MS100S

## SPECIFICATIONS

<b>Model</b>	IV-MS100S
<b>Function</b>	Stirring only
<b>Temp. Range &amp; Accuracy</b>	NA
<b>Heating Power (W)</b>	NA
<b>Plate Type</b>	Acrylic
<b>Display Indication</b>	Time & Speed
<b>Speed Range</b>	200-1500rpm
<b>Maximum Stirring Volume</b>	2000ml
<b>Timing range</b>	0-999min
<b>Display</b>	LCD
<b>Control method</b>	Knob type
<b>Plate size</b>	130mm Ø
<b>Power Supply</b>	AC 220V / 50-60Hz
<b>Dimension &amp; Weight</b>	345 x 225 x160mm & Weight 2kg