

IV-MS5LH Ceramic Hotplate Magnetic Stirrer

- Digital LCD Display.
- Last Measurement in Memory.
- Heating & Stirring can proceed Simultaneously.
- Bushless Motor for Low Noise.
- Anti scratch and Corrosion Resistance.
- Ceramic Plate Surface.
- ABS Cabinet.



SPECIFICATIONS

Speed range	100-2000rpm
Maximum mixing capacity	5000 ml
Temp. Range & Accuracy	300oC Accuracy +1oC
Temperature Sensor	PT100
Time range	0-99h.99min
Top Plat Size	168 x 168mm
Display	LCD
Motor Type	Bushless Motor
Operating Mode	Temperature, Speed, Timing, Memory (last data)
Power	500W
Power Supply	AC 220V / 50-60Hz
Packing Dimension (LxWxH)	335 x 235 x170mm
Net / Gross Weight	2.5kg / 3kg