

IV-OHS100 Overhead Stirrer

- Digital LCD Display.
- 100 ~ 2500rpm speed.
- Maximum Volume 20L.
- Stirring has three modes - Clockwise, Anticlock wise & Alternating.
- Stainless Steel Head.
- Digital type bushless motor.



SPECIFICATIONS

Max. Volume	20L
Display	LCD
Speed(r.p.m)	100~2500
Working Time	Continuous
Structure Exterior	Plastic shell
Motor power	60W
Setting & Speed mode	Digital
Timer	Yes
Max Viscosity	10000 Mpas
Plate size	200 × 300 mm
Plate load	50kg
Power Supply	Using Adaptor, AC 220V / 50-60Hz
Dimension & Weight	400 × 250 × 210 mm & 7.5kg