

## IV- IF100D Ice Flaker

- RED LED Display.
- Fully Automatic microprocessor control.
- Stainless Steel construction.
- Continuous Ice Flakes Output.
- CFC free compressor.
- Low water level Indication.
- Overload protection.
- Noiseless operation.
- Front & Back Ventilated Design.



## SPECIFICATIONS

<b>MODEL</b>	IF100D
<b>Capacity of Ice flake (kg per 24 hrs.)</b>	100 kg
<b>Storage Capacity (kg ~ app.)</b>	25 kg
<b>Display</b>	RED LED
<b>Storage bin</b>	ABS Plastic Bin
<b>Ambient temperature</b>	Below 25 °C for optimum working
<b>Construction</b>	Stainless Steel
<b>Noise level</b>	≤52 dB
<b>Water Supply</b>	3 / 4 inch BSP
<b>Cooling Method &amp; Refrigerant</b>	Air & Refrigerant R290A
<b>Power supply</b>	AC 220V / 50Hz
<b>Input Power</b>	420W
<b>Packing Dimension LxWxH</b>	548X611X883mm
<b>Gross Weight</b>	64kg