

IV-DWB4 Digital Water Bath



- Digital Thermostatic Display.
- High Precise Temperature Control.

SPECIFICATIONS

| | |
|----------------------------|-----------------------------|
| MODEL | DWB4 |
| Column | 2 Rows 4 Holes |
| Internal volume (approx.) | 12L |
| Temp. Range | RT+5~99.9oC |
| Temp. Resolution, Accuracy | 0.1 C, Accuracy ± 0.5 C |
| Water sink | Stainless steel |
| Temp. control mode | PID |
| Cover aperture | 120mm |
| Cover Ring per hole | 5pcs. |
| Power Supply | AC 220V/ 3.6A(50/60Hz) |
| Power rating | 800W |
| Inner Chamber size (LxWxH) | 300 x 300 x 140 mm |
| Packing Dimension (LxWxH) | 410 x 410 x 260 mm |
| Gross Weight | 6kg |