

IV-DWB6 Digital Water Bath



- Stainless Steel Tank.
- High Corrosion Resistance.

SPECIFICATIONS

MODEL	IV-DWB6
Column	2 Rows 6 Holes
Internal volume (approx.)	18L
Temp. Range	RT+5~99.9oC
Temp. Resolution, Accuracy	0.1 C, Accuracy ± 0.5 C
Water sink	Stainless steel
Temp. control mode	PID
Cover aperture	120mm
Cover Ring per hole	5pcs.
Power Supply	AC 220V/ 3.6A(50/60Hz)
Power rating	1200W
Inner Chamber size (LxWxH)	450 x 300 x 130 mm
Packing Dimension (LxWxH)	660 x 400 x 260 mm
Gross Weight	8kg