

IV-DV62T Viscometer - Touch screen

- Touch screen Display.
- Memory upto 30 sets.
- Temperature Monitoring.
- Automatic Detection.
- Selectable Speed.
- Date & Time.
- Product Name Input.
- Optional Accessory - Low Viscosity Adaptor(0 Spindle).



SPECIFICATIONS

Display	Touch screen
Measuring range	50-20,00,000 mPas
Measuring Accuracy	±2 % (Newtonian liquid)
No of speed	8 Speed Option (0.3, 0.6, 1.5, 3, 6, 12, 30, 60 rpm)
Ambient Temperature	10 C to 50 C
No of Spindle(rotor)	4
Zero (0) Spindle(rotor)	Optional 0 Spindle (or 1 to 50mPas)
Interface	RS232C
Memory	30 sets
Date & Time	Yes
Single Measurement Auto Stop	Afer measurement result is stable, spindle stop automatically
Power supply	AC 220V/50Hz
Dimension (LxWxH) & Weight	475 × 200 × 375 mm & 6.8 kg