

IV-PS650 Probe Sonicator - Touch screen

- Touch screen display - 4.3 inches TFT.
- Touch screen display - 4.3 inches TFT.
- User Password Protection available.
- Microprocessor controlled and completely programmable.
- Large LCD display with Time setting .
- Temperature Indicator and controller.
- Auto-tuning for convenience and optimal processing efficiency.
- Variable Power Output control.
- Auto frequency tracking.
- Output amplitude 0-100% adjustment.
- Automatic fault alarm.



SPECIFICATIONS

Frequency	20-25KHz - Automatic tracking
Process capacity	0.5-500ml (approx.)
Standard Probd (Titanium Alloy)	6mm (1/4") Standard
Display	4.3 inches TFT
Display Function	Temperature, Power, Time, etc.
Temperature range(oC)	0-99oC
Timer (Selectable)	Yes
Overload Protection	Yes
User Password Protection	Available
Input method & Data Memory	Touch screen control & 20 groups
Output Power	6.5-650W (adjustable)
Power supply	AC 220V /50HZ
Dimension LxWxH & Weight	490 x 370 x 280mm & 10 kg + 5kg (Sound Proof Enclosure)