

IV-UC SERIES UC3 Ultrasonic Cleaner

- Inner Tank & Outer body made of Stainless Steel.
- Digital Control.
- RED LED display.
- Heating Temperature RT+5 ~ 80°C.
- Temperature with timer function.
- Drainage and cool fan available in 6L to 24L.
- LID & Basket as Standard Supply.



SPECIFICATIONS

MODEL	UC3	
Tank Capacity (App)	2.5L	
Tank Size (mm)	240x140x100	
Overall Size (mm)	270x170x240	
Ultrasonic Power	100W	
Heating power	100W	
Tank Drain	NA	Drain Valve
Time Setting	1-99 min.	
Heating Temp.	RT. +5 ~ 80 C	
Display	Digital RED LED	
Frequency	40kHz	
Inner Tank Material	SS 304	
Outer Tank Material	SS 201	
Weight (approx.)	4 kg	

* Specifications are subject to change without Notice