

## IV-UC SERIES UC12 Ultrasonic Cleaner

- Inner Tank & Outer body made of Stainless Steel.
- Digital Control.
- RED LED display.
- Heating Temperature RT+5 ~ 80°C.
- Temperature with timer function.
- Drainage and cool fan available in 6L to 24L.
- LID & Basket as Standard Supply.



## SPECIFICATIONS

<b>MODEL</b>	UC12	
<b>Tank Capacity (App)</b>	11L	
<b>Tank Size (mm)</b>	330x300x150	
<b>Overall Size (mm)</b>	360x330x310	
<b>Ultrasonic Power</b>	300W	
<b>Heating power</b>	400W	
<b>Tank Drain</b>	NA	Drain Valve
<b>Time Setting</b>	1-99 min.	
<b>Heating Temp.</b>	RT. +5 ~ 80 C	
<b>Display</b>	Digital RED LED	
<b>Frequency</b>	40kHz	
<b>Inner Tank Material</b>	SS 304	
<b>Outer Tank Material</b>	SS 201	
<b>Weight (approx.)</b>	12 kg	

Note : Ultrasonic Cleaner is not recommended for the following gemstones Opal, Pearl, Emerald, Tanzanite, Malachite, Turquoise, Lapis and Coral