

Bottle Top Dispenser - Solvent Resistance

- 5-30ml is suitable for continuously dispensing liquid.
- Solvent Resistance: General Acid and Alkali use.
- Material - PP (Poly Propylene).
- 5 Standard Adaptor dia.

5 Standard Adaptor dia.

Adaptor dia 45 / 40 mm.

Adaptor dia 45 / 32 mm.

Adaptor dia 45 / 38 mm.

Adaptor dia 32 / 25 mm.

Adaptor dia 32 / 28 mm.



SPECIFICATIONS

| | |
|---------------------------|-------------------------|
| Volume | 5-30ml |
| Test Volume | 5ml |
| Accuracy | ± 0.7ml |
| Solvent Resistance | General Acid and Alkali |
| Material | PP |
| Filing Tube | 220mm |
| Bottle | NA |
| Standard Adaptor dia | 5 sizes |
| Density max. | 2.2g/cm ² |
| Adjustable knob control | Yes |
| Packing Dimension (LxWxH) | 290 × 210 × 95 mm |
| Gross Weight | 1.5kg |