



## IV- LTOC1000 Total Organic Carbon Analyzer (TOC)

### SPECIFICATIONS

Make	LABMAN
Detection range	0.001-1.500mg/l (ppm)
Detection limit	1ppb
Max tolerance	±5%
Repeatability tolerance	±3%
Analysis time	< 3 min
Measurement Range	1µ g/L -1000 µ g/L
Response time	< 10 min
Sample temperature	1-95°C
Environment temperature	10-40°C with temperature change ±5 °C/day
Mode	Offline
Software	In-Built
Auto Sampler	Optional
Interface	USB
Relative Humidity	≤ 85%
Display	7" color Touchscreen
Memory	8GB capacity, and no limitation for storage time
Power	100W
Power Supply	AC 220V,50/60Hz
Packing Dimension (LxWxH)	440 x 220 x 300mm