

## IV-128 LABORATORY STIRRER / MECHANICAL AGITATOR

Heavy & Stable cast iron base duly painted with mild steel chrome plated upright rod. AC/DC motor with stainless steel stirring shaft and blades, variable speed control unit provided. 1/35 H.P. motor provided for smooth longtime stirring. Height of stirrer and stirring speed adjustable. To work on 220V AC 50Hz single phase.



## SPECIFICATIONS

<b>Power</b>	220 Volts
<b>Phase</b>	Single Phase
<b>Motor</b>	1/35 HP
<b>Shaft MOC</b>	Stainless Steel
<b>Speed</b>	Adjustable
<b>Base</b>	Heavy Cast Iron