

IV-73 DIGITAL DISINTEGRATION TEST APPARATUS

PRODUCT FEATURES

- Programmable Time With Digital Display
- Auto Stop After Completion of Set Time With Alarm
- Memory for Protection of Running Test in Case of Power Failure
- Separate Regulator for Control of Heating
- Facility for Hold and Resume the Running Test for inspection of Sample

ACCESSORIES

- External Temperature Sensor
- Mains Lead
- Instruction Manual
- 1000 ml Glass Beaker 2 Nos.
- Disc 12 Nos.
- Tube 12 Nos.



Technical specification

Type of Heating	Hot Plate
Up-down Speed	30 ± 2 strokes/min
Up-down Length	55 ± 5mm
Time Control	Digital (Programmable from 1 sec to 99hr 59min, 59 sec)
Electric Motor	Low rpm A.C. Motor
No. of Basket	One
Audio Alarm at completion	Yes
Hold & Resume Facility	Yes
Reaction beaker	1000ml
Resolution	±2%
Power	230V ±10% AC 50hz
Dimension	365x350x495mm (approx)
Weight	9kg (approx)