

IV-71 DIGITAL FRIABILITY TEST APPARATUS

PRODUCT FEATURES

- Programmable Time With Digital Display
- Auto Stop After Completion of Set Time With Alarm
- Auto Removal of Sample Without Opening the Drum
- Precise Control of RPM 25±1

ACCESSORIES

- Acrylic Drum 2Nos
- Main lead
- Instruction manual



Technical specification

LT-71

Model	IV-71
Rotation Speed	25±1
Time Range	1 sec to 99 min 59 sec
No Of Wheels	TWO
Test Time Programmable	Programmable
Electric Motor	Low rpm A.C. Motor
Auto Removal of sample	Yes
Audio Alarm at completion	Yes
Power	230V±10% AC, 50hz
Dimension	270x280x290mm (approx)
Weight	10kg