

## IV-376 MICROPROCESSOR COLONY COUNTERS

### PRODUCT FEATURES

- Count Correction Facility
- Result-Count Storage
- Audible Confirmation of Each Count
- Auto Marker Pen
- Uniform Glare Free Illumination
- Wolffhugel Glass Grid with Focussing Facility
- LED Based Illumination

### ACCESSORIES

- Counting Pen
- Magnifier Glass
- Stand
- Instruction Manual



### Technical specification

<b>Display</b>	6- Digit Seven Segment LED
<b>Max Count</b>	999999
<b>Keyboard</b>	6 soft touch membrane keys
<b>Correction</b>	Count correction facility available to cancel the unwanted counts
<b>Average</b>	Average mode facility available
<b>Dish Size</b>	110 mm
<b>Magnification</b>	x 1.7 surface magnification
<b>Reset</b>	Data reset facility available
<b>Memory</b>	99 memory location
<b>Storage</b>	99 result count storage retained even if the system is switched off
<b>Dimensions (L x B x H)</b>	255 x 220 x 100 mm (Approx.)
<b>Weight</b>	3 Kg. (Approx.)
<b>Power</b>	230 V±10% AC, 50 Hz