

IV-37 MICROPROCESSOR COLONY COUNTERS

PRODUCT FEATURES

- Count Correction Facility
- Result-Count Storage
- Audible Confirmation of Each Count
- Auto Marker Pen
- Uniform Glare Free Illumination
- Wolffhugel Glass Grid with Focussing Facility
- LED Based Illumination

ACCESSORIES

- Counting Pen
- Magnifier Glass
- Stand
- Instruction Manual



Technical specification

Display	4- Digit Seven Segment LED
Max Count	9999
Keyboard	6 soft touch membrane keys
Correction	Count correction facility available to cancel the unwanted counts
Storage	40 Result-count storage retained even if the system is switched off
Average	Average mode facility available
Dish Size	110 mm
Magnification	x 1.7 surface magnification
Dimensions (I X B X H)	255 x 220 x 100 mm (Approx.)
Weight	3 Kg, (Approx.)
Power	230 V±10% AC, 50 Hz