

IV-21 DIFFERENTIAL CELL COUNTER

PRODUCT FEATURES

- Facility to count 8 type of cells
- 16x2 line alphanumeric LCD display
- Automatic percentage (%) calculation
- Individual cell counting facility
- Audible beep facility

ACCESSORIES

- Dust cover
- Instruction manual



Technical specification

Cell Counting Key	16 keys
Total Keys	16 keys (8 cell keys and 3 function keys)
Body	Metal body
Display	16x2 line LCD display
Mains Input	230±10% ac, 50 hz
Weight	1.5 kg
Dimensions	250 x 230 x 150 mm (L X B X H)