

IV-63 AUTO KARL FISCHER TITRIMETER

PRODUCT FEATURES

- Auto Titration Control
- Auto Zero Burette
- Adjustable Timer
- Visual & Audible Alarm

ACCESSORIES

- Dual platinum electrode
- Burette clear glass 10ml
- Rubber bellow
- Titration Vessel
- Reagent Bottel



Technical specification

Display	Array of 10 LEDS
Sample Size	1-50 mg of H ₂ O
Indication	RUN, ADD, END LEDS
Kf Flow	Auto Control Valve
Burette	Auto Zero 10ml Burette
Timer	Adjustable Timer
Dispensing Resolution	0.05ml
Measuring Electrode	Dual Platinum KF Electrode
Power	230 V \pm 10% , 50 hz , AC