

IV-430 DIGITAL BLOOD ROLLER MIXER

PRODUCT FEATURES

- Rocking and rolling action for complete mixing
- six roller design with small footprint
- designed for continous quiet operation
- adjustable speed range from 10 to 70 rpm



ACCESSORIES

- Main lead
- Dust cover
- Instruction manual

Technical specification

Display	LCD
Adjustable speed range	10 to 70 rpm
No. of Roller	6
Maximum load	10 test tube
Timer	Yes
Speed	Adustable
Weight	3kg (approx)
Power supply	230V±10% ac 50Hz