

IV- 44 ROTARY BLOOD MIXER

PRODUCT FEATURES

- AC Synchronous motor for noiseless operation.
- Suitable for a variety of centrifugal tube 1.5-50ml.
- Circular rotating action gives gentle and effective mixing.
- Fix speed of 30 RPM with an accuracy of ± 2 RPM.
- Continuous/infinite mode of operation.
- Fix plate angle of 45 Degree.
- Easy to operate, long running, stable and low noise.
- Maximum operational weight capacity of 1.5 kgs.
- Easy to remove and install disk plates for replacement and cleaning.



Technical specification

Speed	30 RPM
Slope angle	45°
Voltage	220 V AC
Power	6V
Maximum operational weight capacity	1.5 kg
Dimension(WDXH) MM	165X150X145