

IV- 780 DIGITAL TABLET HARDNESS TESTER



Technical specification

Display	LCD with backill display
Keyboard	User Friendly 9 Keys operations
Hardness Range	5-500 N
Hardness Accuracy	± 0.05% FS
Measuring Units	N(Newton), KGF (Kilogram Force) and Lbf (Low bound Force)
Interface Ports	RS-232
Battery	1200 mAH, 8.4V, Ni-Hi Rechargeable Battery
Peak Force Measurement	Yes
Acrylic Safety Glass Cover	Yes
Stored Results	Upto 10 (also calculate the average value of test results)
Tablet Size	Yes (2mm to 28 mm)