

## IV-113 DIGITAL HAEMOGLOBIN METER

### ACCESSORIES

- Glass test tubes - 5 Nos.
- Instruments cover - 1 Nos.
- Operating Manual - 1 Nos.
- Spare Lamp - 1 Nos. (Fitted inside)



### Technical specification

<b>Display</b>	3 Digit Red LED Display
<b>Display Modes</b>	%T.O.D. and 3 Concentration Modes C1, C2, C3
<b>Zero Adjustment</b>	Automatic (soft touch key)
<b>Accuracy</b>	± 0.02 O.D.
<b>Stability</b>	± 0.02 O.D. Per Hour
<b>Light Source</b>	3.5V LED bulb
<b>Amplifier</b>	Micro Controller Based.
<b>Memory</b>	In built memory to retain last "K" factor for concentration even after power loss.
<b>Sample Volume</b>	1 ml. Minimum.