

IV-110 DIGITAL DRY MELTING POINT APPARATUS


PRODUCT FEATURES

- Furnace based dry heating
- Highly accurate and stable
- Excellent view angle
- Fast heating

ACCESSORIES

- Operation Manual
- Capillary Tubes
- Power Cable
- Test Tube for Boiling Samples

VISUAL-WITHOUT SILICON OIL



Technical specification

Display	3 Digit LED
Display Format	Menu driven
Setting	By Soft Touch Keys
Heating Devices	Dry Furnace
Temp. Sensing Devices	RTD/PT100
Range	5°C above ambient to 350°C
Resolution	1°C
Accuracy	Depend on manual viewing/operation
Calibration	Manual
Power Supply	230 V \pm 10% AC 50 HZ
Dimension	355x340x450mm (approx)
Weight	6kg