

IV-109 MICROPROCESSOR MELTING POINT APPARATUS

PRODUCT FEATURES

- 128x64 Graphical LCD display with backlight
- Icon-based navigation menu for ease of operation
- Selectable heating modes (slow, medium and Rapid)
- Storage up to 1000 results
- Quick Cooling with Integrated fans
- Audio and Visual Indicators for System Status
- Quick Cooling with integrated Fans

ACCESSORIES

- Capillaries 1 Box
- Cleaning tool 1 Nos
- USB Cable 1 Nos
- Software 1 Nos

AUTOMATIC WITHOUT SILICON OIL



Technical specification

No of sample	1 sample at a time
Temperature range	Ambient to 300°C
Temperature accuracy	±1°C at 200°C and ±2.5°C at 275°C
Temperature sensor	Platinum resistance PT-100
Selectable heating mode	Slow, medium and rapid
Display	128x64 graphical LCD display
No. of key	4
Memory	Storage up to 1000 results
Power Supply	230 V ±10% AC 50 HZ
Dimension	300x200x110mm (approx)
Weight	3kg