

IV-77 BULK DENSITY APPARATUS

PRODUCT FEATURES

- Extremely useful for Pharmaceutical Lab
- Highly Accurate and Easy to Operate
- Microcontroller Based Digital Display
- Soft Touch Membrane Keys
- 3 LED's for status Indication
- Audible Alarm Sound
- Last Count Retention on Power Failure

ACCESSORIES

- Glass cylinder 2 Nos (100ml or 250ml)
- Cylinder holder
- Main lead
- Instruction manual



Technical specification

Maximum Counts	9999
Display	4 Digit 7-Segment Red LED
Keyboard	6 Soft Touch Membrane Keys
Cylinder	2 Nos., 100ml/250 ml Cylinders
Capacity of cylinder	100ml/250ml (mention during orders)
Drop Height	25 mm (Approx.)
Power	230 V \pm 10% AC, 50 Hz.
Dimensions	600 x 450 x 200 mm (Approx.)
Weight	12 Kg (approx.)