

## IV-175 MICROPROCESSOR AUTO LEAK TEST APPARATUS

### PRODUCT FEATURES

- Programmable Time With Digital Display
- Auto Step After Completion Of Time With Alarm
- Precise Control Of Vacuum Pump
- Heavy Duty Oil-Free Vacuum Pump
- 1 mm Thick SS-304 Cabinet

### ACCESSORIES

- Instruction manual
- 12 Inch vacuum desiccator
- Connecting pipes

12" DESICCATOR



### Technical specification

<b>Type</b>	GMP Model, In-built Test Unit
<b>Time Control &amp; Range</b>	Digital, Programmable, 0.01 sec to 99 min 59 sec
<b>Test Vacuum Range</b>	Adjustable up-to 650 mm Hg
<b>Vacuum Control Range</b>	Digital Display
<b>Vacuum control Release</b>	By Single Front Panel Knob
<b>Audio Alarm at Completion</b>	Yes
<b>Vacuum Pump</b>	Oil free, MV 27" $\pm$ 5%, 1/4HP Motor
<b>Test Unit Size</b>	200mm, 250mm & 300mm (any one) (Specify in Order)
<b>Power</b>	230V $\pm$ 10% AC 50Hz
<b>Dimension</b>	600x450x200mm (approx)
<b>Weight</b>	8kg (approx)