

IV-160 DIGITAL BOMB CALORIMETER

PRODUCT FEATURES

- Manual measurement & calculation.
- Initial temperature rise & fall display.
- Accurate result

ACCESSORIES

- Steel drum
- Bomb with sample holder
- Stirrer
- Temperature Probe
- Nichrome Wire
- Pallet press
- Cotton thread
- Safely valve for high pressure
- Oxygen pressure gauge
- Benzoic acid powder
- Tool kit
- Operation manual
- Power supply cable (Main lead)



Technical specification

Display	3 Digit LED
Temperature accuracy	0.1°C
Temperature Resolution	0.1°C
Temperature Detector	RTD
Temp. Range	RT. to 99.9°C
Calculate	Manual
Power supply	230V±10%/50Hz