

IV-350 ULTRASONIC CLEANERS

PRODUCT FEATURES

- Stainless steel tank has resistance to wear and long work life.
- With digital timer & heater control.
- With Stainless steel basket
- Use just tap water, or industrial alcohol and solvent cleaner for higher cleaning requirement
- Industrial control chip micro controller flexible circuit boards control
- Secure & stable
- Temperature tunable from 20°C to 80°C.

ACCESSORIES

- SS Basket
- Main Lead
- Instruction Manual



Technical specification

Ultrasonic Frequency	40,000 Hz
Material	Stainless Steel SS304
Tank Capacity	2.0L, 3.5L, 4.5L, 6L
Timer	1-30 minutes, Digital timer
Power Supply	AC 100-120, 50/60Hz AC 220-240 50/60 Hz 240W
Ultrasonic Power	240W
Heating Power	200 W, Digital Heater
Unit Size	325 x 265 x 280 mm
Tank Size	300 x 240 x 150 mm
Package Size	430 x 365 x 420 mm
Net. Weight	6.0 kg
Gross. Weight	7.0 kg