

IV-101 DIGITAL PH METER

PRODUCT FEATURES

- Inbuilt stirrer
- Bright LED display
- Stable and accurate result
- Auto temperature compensation

ACCESSORIES

- Combination pH Electrode
- Buffer Tablets 4 pH & 7 pH
- Operation Manual
- Buffer Bottles (4 pH & 7 pH)
- Dust Cover
- Temperature Probe
- Rotational pH Electrode Stand

INBUILT STIRRER



Technical specification

Parameter	pH	mv
Range	0.00 to 14.00 pH	0 to 1999.9 mV
Resolution	0.01pH	1 mV
Accuracy	±0.01 pH	±1mV
Display	3 1/2 Digit LED	
Temp Compensation	Manual: 0 to 100	
Input Impedance	> 10 ohms	
Slope Control	80 to 120%	
Recorder Output	0 to 10 mV/pH Adjustable 0 to 10 mV/100 mV Adjustable	
Power	230 V ±10% AC, 50Hz	
Dimensions	275 x175x 76 ctrn (LxBxH)	
Sensor	Combined pH Electrode	
Stirrer	inbuilt Stirrer (with speed variable)	