

## IV-23 DIGITAL PH-CONDUCTIVITY TEMPERATURE METER

### PRODUCT FEATURES

- 3 parameters in single instrument (pH/Conductivity/temperature)
- Compact design
- Bright LED display
- Accurate and Stable

### ACCESSORIES

- Conductivity cell 1.0
- Combination pH Electrode
- Buffer Tablets 4 pH & 7 pH
- Operation Manual
- Dust Cover
- Rotational pH Electrode Stand
- Temperature probe



### Technical specification

#### pH Measurements

<b>pH</b>	0 to 14 pH
<b>Resolution</b>	0.01 pH
<b>Accuracy</b>	$\pm 0.01$ pH $\pm 1$ digit

#### Conductivity Meter

<b>Range</b>	0 to 2 & 20mMhos/cm
<b>Accuracy</b>	0.5% $\pm 1$ digit

#### Temperature Measurement

<b>Range</b>	0 to 100°C
<b>Resolution</b>	0.1°C
<b>Accuracy</b>	0.2% + 1 digit
<b>Display</b>	3½ digit LED display
<b>Power</b>	230V $\pm 10\%$ AC 50hz
<b>Dimensation</b>	275x125x175mm approx
<b>Weight</b>	3kg